

Products Promotion: Portable & Lab Meter

Pocket Tester with replaceable probe, Waterproof pH ~ °C ~ ORP ~ Conductivity ~ TDS



pH55 pH/Temp Meter, Range: -2.0 to 16.0pH / -5.0 to 60.0°C; Accuracy: ±0.1pH, ±0.5°C

ORP57 ORP/Temp Meter, Range: ±1000mV / -5.0 to 60.0°C / Accuracy: ±2mV, ±0.5°C

EC59 Conductivity/TDS/Temp Meter, Range: 3999µS/cm / 2000ppm / 0.0 to 60.0°C

EC60 Conductivity/TDS/Temp Meter, Range: 20.00mS/cm / 10.00ppt / 0.0 to 60.0°C

*All meters are supplied complete with protective cap, carton box, batteries and instruction

Economic Portable Meter pH ~ °C ~ ORP ~ Conductivity ~ Resistivity ~ TDS

MW101 pH Meter, Range: 0.00 to 14.00pH; Accuracy: ±0.02pH

MW102-1 pH/Temp Meter, Range: -2.00 to 16.00pH / -5 to 70°C; Accuracy: ±0.02pH, ±0.5°C
-Come with double junction standard pH electrode & temp. probe

MW500 ORP Meter, Range: ±100mV; Accuracy: ±5mV

SC110 Conductivity/Resistivity/Salinity/TDS/Temp Meter, 450 sets memory
-Range: 0.00µS.cm⁻¹ to 200.0mS.cm⁻¹ / 0.00 to 100.00MΩ.cm / 0 to 70ppt / 0 to 2000ppm

MW600 Low Cost DO Meter, Range: 0.0 to 19.9 mg/L, Accuracy: ±1.5%FS

Mi605 DO meter, Range: 0.0 to 45.00ppm, Accuracy: ±1.5%FS

*All meters are supplied complete with electrode, solutions, batteries and manual.



Digital Refractometer Brix ~ Salinity

MA871 Brix Meter, Range: 0 to 85% Brix / 0 to 80°C; Accuracy: ±0.2% , ±0.3°C

Refractometer for Fructose, Glucose, Invert Sugar. Range: 0 to 85% mass, Accuracy: ±0.2%

MA872
Fructose

MA873
Glucose

MA881
Invert Sugar

* Supplied in a carton box, complete with 9V battery and manual.



Products Promotion: Portable & Lab Meter

Chlorine Photometer

Free ~ Total Chlorine



MW10 Budget Free Chlorine Meter, Range: 0.00 to 2.50 mg/L

MW11 Budget Total Chlorine Meter, Range: 0.00 to 3.50mg/L

Mi404 Portable Free & Total Chlorine Meter
- Range: 0.00 to 5.00mg/L

*MW10& MW11: Include 2 cuvettes, reagents and manual

*MW404: Include 2 cuvettes, powder reagents for FCI & TCI (100 sachets each)



Pocket-Sized Thermometer, Accuracy: $\pm 0.5^{\circ}\text{C}$

$^{\circ}\text{C}$



TH310 Stick Thermometer, Range: -50.0 to 150.0°C

TH300 Portable Thermometer, Range: -50.0 to 150.0°C
-Probe comes with 1m cable.



Laboratory Meter

pH ~ Conductivity ~ DO ~ Water Still ~ Hotplate Stirrer

SP-2100 pH/ORP/Temp Meter, Range: $-2.00 \sim 16.00\text{pH}$ / $-10.0 \sim 110.0^{\circ}\text{C}$
*Include glass pH probe, temp. probe, electrode holder, calibration buffer.

Mi180 pH/mV/EC/TDS/NaCl/Temp. Bench Meter
*Include pH probe, temp. probe, cond. Probe, electrode holder, software, calibration solution.

Mi190 DO/Temp Bench Meter, Range: 0 to 45ppm, -5 to 55°C
*Include DO probe, electrolyte solution, software.

CC-570 Colony Counter, 4 digit readout, for 10cm to 15cm dish

SH-301 Hotplate Stirrer, Range: Up to 450°C , 60-1100rpm, with safety protection

WS400 Auto Water Still with conductivity display, 4L/hour, with safety protection
*Also available in 4L/hour double distilled & 8L/hour models.

