

Hotplate & Stirrer

SH-301

Features:

- ◆ Easy and user-friendly operation.
- ◆ Precise heating and accurate temperature control.
- ◆ All in one glass-ceramic top plate: easy clean, chemical-resistant, corrosion-resistant, and scrape-resistant.
- ◆ LED Indicators for Power, Heating and Stirring.
- ◆ Safety first design and performance:
 - 1) Hot warning indicating light: Automatic warning when the temperature of the ceramic top is over 55°C.
 - 2) Overheated safety device: Automatic heat cutoff when the heating temp is over 550°C.
 - 3) Splash resistant design: Avoid sample permeating from the plate or switch.
- ◆ Special circuit board design for power saving and efficient power distributing.
- ◆ Microprocessor adjustable speed control from 60 - 1100 rpm.
- ◆ Advanced motor for extreme smooth rotation, low noise and long life span.



Hotplate & Stirrer

SH-301



Specifications:

Model		SH-301	SH-301A
Function		Heating / Stirrer	
Top Plate	Material	Glass Ceramic	
	Dimension	200 x 200 mm	
Speed	Range	60 - 1100 rpm	
	Control	Adjustable	
Temperature	Range	Ambient temperature to 450°C	
	Control	Adjustable	
Heating Capacity		Max. 750 W	
Hot warning indicator		Yes	
Heating indicator		Yes	
External Temp Probe		No	
Ambient Temperature		5 - 40°C	
Overall Dimension		320 x 210 x 105 mm (W x D x H)	
Weight		3.6 kg	
Power Supply		110VAC ± 10%, 50/60Hz	220VAC ± 10%, 50/60Hz

* The specifications and appearance of the instrument are subject to change without notice.

Agent



SUNTEX INSTRUMENTS CO., LTD.

Web site: www.suntex.com.tw E-mail: suntex@ms1.hinet.net

13F, No. 31, Lane 169, Kang Ning St., Hsi-Chih City, Taipei County, Taiwan, R.O.C. 221

Tel: 886-2-2695-9688 Fax: 886-2-2695-9693

Taichung Branch: 886-4-2319-7996 Kaohsiung Branch: 886-7-334-4103